

SAFILO ANTI COVID-19 PRACTICES

THE GOOD NEWS IS THAT SANITIZING YOUR GLASSES IS EASY



GET PREPARED

Wash your hands frequently - as suggested by healthcare authorities. Wear your face mask and/or shield, latex gloves, and sanitize them after any interaction with the customer and product.



CLEAN ENVIRONMENT

Personally hand each frame to the customer for try on, keeping them separate from the rest of the frames. After the customer visit, restore the set up by sanitizing each item used during the patient engagement including mirrors, shelves and displays.



CLEANING FRAMES PROCEDURES

Sanitize every frame touched using lukewarm water and a gentle pH neutral soap, free of lotion: massage the soap onto each lens, nose pads and earpieces, rinse and dry using a microfiber cloth (do not use paper towels). The use of an ultrasonic bath with this solution is optimal. The cleaning cloth can also be machine washed and used again. In the absence of access to running water, frames can be cleaned with hydrogen peroxide wipes. If hydrogen peroxide wipes are not available, a 50-50 mix of 3% hydrogen peroxide and water solution can be used in a spray bottle to clean the frames with a cleaning cloth. Please be aware that damage to the frame may occur if the solution is not mixed properly in the correct proportions. Remember the best option is the use of soap and water.



DON'Ts

- Do not use rubbing alcohol to sanitize your glasses
- Avoid using products with high concentrations of acid, as they can damage the product (i.e NO acetone, Don't consider hydrogen peroxide as main cleaning agent)
- Do not use other not specific sanitizing wipes, as they can damage the product

Safilo®

THESE INSTRUCTIONS FOR SANITIZATION ARE PROVIDED BY SAFILO AND MEANT FOR ITS EYEWEAR PRODUCT ONLY. The information in this document shall not constitute medical advice. Please always follow the applicable laws and regulations aimed at containing the spread of Covid-19 and consult the local healthcare authorities communications for the most recent developments and advice.